


Redacted

Purpose: To prepare gels for 49143 and 49143 NHA

Procedure: Same as described on page 115.

Samples were resolved on 12% PAGE and then silver stained.



Conclusion 49143 NHA show an enrichment for the 19 kDa protein, and an additional protein below 18.5 kDa (A known protein is about 16 kDa) suggesting we may have selected for mucinase protease. Interestingly 49143 NHA lacks the >106 kDa protein that we initially considered as a potential hemagglutinin. This is very suggestive that the >106 kDa protein could actually be the hemagglutinin!

Redacted

(KPT)